LSCI 3000 Synthetic Biology Workshop Lab Log 3

Chan Cheuk Ka (1155174356) 10-03-2023

1. Experiment 1

1.1. Objective

• To prepare a welled agarose gel for gel electrophoresis.

1.2. Precautions

- Wear gloves.
- Handle the agarose gel with care to prevent breakage.

1.3. Materials

- 0.16g of agarose gel powder
- 20mL of TAE buffer (1×)
- 1μL of RedSafe DNA staining dye (20000×)

1.4. Equipment

1.4.1.Machineries

- Electronic balance
- Microwave

1.4.2. Apparatus

- Spatula
- Pipette (P10) and appropriate pipette tips

1.4.3. Containers

- Conical flask
- Gel casting tank (from the gel maker set)
- Gel tray (from the gel maker set)
- Gel comb (from the gel maker set)

1.5. Preparations

• Set up the gel maker set. Mount the gel tray into the tank.

1.6. Procedures

- 1. Weigh and add 0.16g of agarose gel powder into the conical flask.
- 2. Pour 20mL of TAE buffer into the flask.
- 3. Heat the solution in a microwave for 30 seconds or until all the powder is completely dissolved.
- 4. Cool the flask down to ~60°C in the sink with running water.
- 5. Pipette 1µL of DNA staining dye into the flask.
- 6. Gently mix the sample by swirling.
- 7. Pour the mixture onto the gel tray. Remove air bubbles with a pipette tip if any are formed.
- 8. Insert the gel comb to create $25\mu L$ wells in the gel.
- 9. Wait 30 minutes or until the gel completely solidifies.
- 10. Remove the comb carefully.

2. Experiment 2 (Prerequisite: Lab 2 Experiment 2)

2.1. Objective

• To extract plasmids from the transformed E. coli (DH5α) small cultures.

2.2. Precautions

- Wear gloves.
- Dispose of the pipette tips after each transfer.
- When placing tubes into the centrifuge, position the tubes such that their hinges face outwards to allow consistent pellet settling behaviour.

2.3. Materials

2.3.1.<u>Small cultures</u>

- 2 pRSFDuet-1 transformed DH5α E. coli small cultures in LB broth in falcon tubes
- 2 pLadder6K transformed DH5α E. coli small cultures in LB broth in falcon tubes

2.3.2.iNtRON easy-spin total RNA extraction kit

- 1000µL of resuspension buffer
- 1000µL of lysis buffer
- 1400µL of neutralisation buffer
- 2800µL of washing buffer B

2.3.3. Miscellaneous

- 200µL of ddH₂O
- Ice

2.4. Equipment

2.4.1. Machineries

- Centrifuge
- -20°C refrigerator

2.4.2. Apparatus

• Pipettes (P200/P1000) and appropriate pipette tips

2.4.3. Containers

- 8 centrifuge tubes
- 4 spin columns (from the iNtRON kit)
- 4 collection tubes (from the iNtRON kit)
- Ice box

2.4.4. Miscellaneous

• Marker pen

2.5. Preparations

- Add RNase A solution to the resuspension buffer as per the instructions on the bottle, then chill the buffer with ice until use.
- Add ethanol to washing buffer B as per the instructions on the bottle.
- Centrifuge the E. coli small culture before use.
- Label all centrifuge tubes and columns appropriately.
- Mount the columns into the collection tubes as necessary before use.

2.6. Procedures

- 1. Pipette 250µL of resuspension buffer to each small culture falcon tube.
- 2. Gently mix the samples by vortexing until all bacteria clumps are resuspended.
- 3. Transfer all samples into separate centrifuge tubes with a pipette.
- 4. Pipette 250μL of lysis buffer to each tube.
- 5. Gently mix the samples by inverting the tubes 10 times.

- 6. Let the samples sit at room temperature for 3 minutes.
- 7. Pipette 350µL of neutralisation buffer to each tube.
- 8. Gently mix the samples by inverting the tubes 10 times.
- 9. Let the samples sit in ice for 5 minutes.
- 10. Centrifuge the samples at 13000rpm for 10 minutes.
- 11. Pipette the supernatant fluids into separate columns. Avoid transferring the pellets.
- 12. Centrifuge the samples at 13000rpm for 1 minute.
- 13. Discard the filtrate fluids in the collection tubes and remount the columns into the same collection tubes.
- 14. Pipette 700µL of washing buffer B to each column.
- 15. Centrifuge the samples at 13000rpm for 1 minute.
- 16. Discard the filtrate fluids in the collection tubes and remount the columns into the same collection tubes.
- 17. Centrifuge the samples at 13000rpm for 1 minute or until the filter membranes are dry.
- 18. Put the columns into separate centrifuge tubes. Remove and discard the used collection tubes
- 19. Pipette 50μL of ddH₂O to each column.
- 20. Let the samples stand for 1 minute.
- 21. Centrifuge the samples at 13000rpm for 1 minute.
- 22. Discard the columns and store the plasmids at -20°C for future use.

3. Experiment 3 (*Prerequisite: Lab 2 Experiment 1 + Lab 3 Experiment 1*)

3.1. Objective

• To isolate the desired DNA (GFP encoding gene) from the PCR products with gel electrophoresis.

3.2. Precautions

- Wear gloves.
- Dispose of the pipette tips after each transfer.
- Handle the agarose gel with care to prevent breakage.
- Do not run gel electrophoresis with >140V or the gel could melt.
- Ensure the electrodes of the gel electrophoresis machine are connected correctly before starting the machine.
- Slightly lift the pipette tip while loading samples into the gel wells to minimise spillage.
- Do not look directly into the UV light without the protective shield.

3.3. Materials

- 4 centrifuge tubes each with 25µL of PCR products (3 positive + 1 negative control)
- 6µL of dyed 1kb Plus DNA ladder
- 10μ L of DNA loading dye ($10\times$)
- 1 agarose gel with 25μL wells in a gel tray
- Ample TAE buffer (1×)

3.4. Equipment

3.4.1. Machineries

- 120V constant voltage supply
- UV transilluminator

3.4.2. Apparatus

• Pipettes (P10/P200) and appropriate pipette tips

3.4.3. Containers

- Gel electrophoresis tank
- Gel tray (from the gel maker set)

3.4.4.Miscellaneous

• Paper towels

3.5. Preparations

- Pre-chill the PCR products in ice.
- Put a plastic wrap on the UV transilluminator sample surface before use.

3.6. Procedures

- 1. Put the gel with its gel tray into the gel electrophoresis tank and orient it such that the wells are closer to the black negative electrode.
- 2. Slowly pour TAE buffer into the tank until the gel is completely submerged.
- 3. Carefully pipette 6µL of DNA ladder into the leftmost well of the gel.
- 4. Pipette 2.5μL of DNA loading dye into each PCR product mixture.
- 5. Gently mix the samples.
- 6. Carefully pipette 25µL of each mixture into separate wells of the gel.
- 7. Connect the electrodes of the tank to the constant voltage supply.
- 8. Electrophorese at 120V for 20 minutes or before the samples run off the gel. Check that bubbles are formed at the negative electrode.
- 9. Remove the gel from the tray and transfer it onto a paper towel to dry slightly.
- 10. Place the gel on a UV transilluminator to observe the results.

3.7. Results

The most prominent bands can be observed at ~750bp, which highly correlates to the desired 730bp GFP encoding gene.




Figure 1 - Electrophoresed gel under UV light From left to right: DNA ladder, null, samples 1-3, negative control

4. Experiment 4 (*Prerequisite: Lab 3 Experiment 3*)

4.1. Objective

• To extract and purify the desired DNA (GFP encoding gene) with gel extraction.

4.2. Precautions

- Wear gloves.
- Dispose of the pipette tips after each transfer.
- Do not look directly into the UV light without the protective shield.

4.3. Materials

4.3.1. Electrophoresed gel

• 1 electrophoresed gel

4.3.2.iNtRON MEGAquick-spin plus total fragment DNA purification kit

- $<600\mu L$ of BNL buffer
- 750μL of washing buffer

4.3.3. Miscellaneous

• $40\mu L$ of ddH_2O

4.4. Equipment

4.4.1.Machineries

- UV transilluminator
- 55°C heat block
- Centrifuge
- -20°C refrigerator

4.4.2.Apparatus

- Pipettes (P200/P1000) and appropriate pipette tips
- Scalpel

4.4.3. Containers

- 2 centrifuge tubes
- 1 spin column (from the iNtRON kit)
- 1 collection tube (from the iNtRON kit)

4.4.4.Miscellaneous

• Marker pen

4.5. Preparations

- Add ethanol to the washing buffer as per the instructions on the bottle.
- Put a plastic wrap on the UV transilluminator sample surface before use.
- Label the centrifuge tubes and columns appropriately.
- Mount the column into the collection tube as necessary before use.

4.6. Procedures

- 1. Place the gel on a UV transilluminator.
- 2. Locate a band of the desired length (730bp) with the most prominence and isolate the entire band with a sterile scalpel.
- 3. Weigh the gel slice. Trim until it is \leq 300mg if necessary. (We had 240mg in practice)
- 4. Transfer the gel slice into a centrifuge tube.
- 5. Pipette (weight in mg \times 2) μ L of BNL buffer to the tube. (We pipetted 420 μ L in practice)
- 6. Gently mix the sample by vortexing.
- 7. Place the tube in a heat block at 55°C for 7 minutes or until the gel slice is completely dissolved. Vortex the tube every 2-3 minutes.
- 8. Allow the sample to cool to room temperature.

- 9. Transfer all of the sample mixture to a column with a pipette.
- 10. Centrifuge the sample at $11000 \times g$ for 30 seconds.
- 11. Discard the filtrate fluid in the collection tube and remount the column into the same collection tube.
- 12. Pipette 750μL of washing buffer to the column.
- 13. Centrifuge the sample at $11000 \times g$ for 30 seconds.
- 14. Discard the filtrate fluid in the collection tube and remount the column into the same collection tube.
- 15. Centrifuge the sample at $16000 \times g$ for 3 minute or until the filter membrane is dry.
- 16. Put the column into a sterile centrifuge tube. Remove and discard the used collection tube.
- 17. Pipette 40μL of ddH₂O directly onto the filter membrane of the column.
- 18. Let the sample stand for 1 minute.
- 19. Centrifuge the sample at $16000 \times g$ for 1 minute.
- 20. Discard the column and store the purified DNA at -20°C for future use.